SCUDO

TOTAL CONTACT DIABETIC CAST

A TOTAL CONTACT INNOVATION

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The technological
revolution
that leans
on science.

FGP E MEDI 5
FGP AND MEDI

TOTAL CONTACT CASTING 6
TOTAL CONTACT CASTING


IL FUTURO. SCUDO, TOTAL CONTACT DIABETIC CAST 12
THE FUTURE. SCUDO, TOTAL CONTACT DIABETIC CAST

TEST DI LABORATORIO 19
LABORATORY TEST

SCHEDA TECNICA 22
DATA SHEET

RICAMBI 24
SPARE PARTS

BIBLIOGRAFIA 26
REFERENCES

FGP

Orthopedic
& Rehabilitation
Equipment

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FGP and MEDI. Two protagonists of well-being with a common vision.

At FGP we believe that freedom of movement is quality of life.

*For this reason, since 1988, **we are leader in Italy in the design, production and distribution of orthopaedic and rehabilitation devices.***

Competence, dynamism and innovation makes FGP a brand capable of anticipating the needs of patients, through dialogue with the scientific world specialists and technicians.

*The acquired know-how in the orthopedic field and the important collaborations **with diabetic foot experts** have led to the development of high-performance products based on solid scientific bases for a so specific treatment area.*

At the beginning of 2023 FGP joined the international group MEDI, one of the leading manufacturers of elastic stockings and health products in the world and already a shareholder of the company: a path of collaboration was created by sharing vision and values, strengthening an already unique know-how in the sector that places FGP among the most innovative companies in the production of orthopaedic braces.

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Diabetic foot ulcer overview

Neuropathic DFUs are chronic wounds that result from a combination of factors. Some of the main factors that result in DFUs include:

- 1 Excessive direct pressure secondary to deformities** (the weight on the feet that naturally results from standing);
- 2 Shear** (usually twisting forces between the skin and footwear);
- 3 Repetitive stresses** (for example, from walking);
- 4 Peripheral neuropathy** (damage to the nerves in the feet).

These factors rarely cause problems in healthy people but the accompanying neuropathy means that the "pain feedback loop" that induces humans to limp or keep their weight off an injury no longer functions properly in diabetic patients. The result is that pressure and abrasion injuries, such as bruising and blisters caused by poor footwear, are not noticed and often ignored, and the continual use of the limb disrupts and damages the healing process.

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Explanation of off-loading techniques

*Effective off-loading is not simple, especially because of the extensive healing time often required. Most patients need to **remain mobile to continue with daily life**, therefore enforced restrained bed rest are neither convenient nor recommended method to achieve effective off-loading. Other off-loading techniques such as crutches and wheelchairs are not supported by published clinical evidence as effective methods.*

OFF-LOADING DEVICES, THAT ARE ACTUALLY USED, CAN BE FUNCTIONALLY CLASSIFIED INTO TWO GROUPS:

REMOVABLE DEVICES, which let the patient control the amount of use;

IRREMOVABLE DEVICES, which are not patient-removable and must be removed or changed by health professionals.

A recent Cochrane review of off-loading devices found that non-removable, pressure-reducing casts were the most effective off-loading devices for DFUs.

Immobilization of the ankle joint reduces pressure in the forefoot and plantar regions during the propulsion gait phase, which occurs during forward motion. Both TCCs and removable cast walkers can immobilize the ankle, but TCC ensures constant use by the patient.

TCC reduces the stride length, which slows down the foot and reduces the force applied to the foot. TCC has been shown to reduce plantar pressure by 32%, 63%, and 69% on the fifth, fourth, and first metatarsal heads, respectively; 65% on the great toe; and 45% on the heel.

TCC can be used to off-load the heel, but a panel of experts agreed that this approach might not be as effective as for off-loading of forefoot ulcers.

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TOTAL CONTACT CASTING (TCC) is a specially designed cast, designed to reduce plantar pressure peaks by redistributing loads over a larger surface (off-loading) in patients with ulcerative lesions (DFUs) secondary to neuroischemic complications of diabetes mellitus. Reducing the pressure on the lesion has been shown to be an essential therapeutic component for healing of meso-retrofoot plantar lesions.

The use of total contact cast (TCC) requires a specific knowledge and manual skills to prepare a structure, usually made in fibreglass, that includes the foot (including the fingers) and the leg (up to the middle third); this device allows to redistribute the pressures from the foot on the leg. This pressure reduction removes the the lesion traumatism component allowing a more physiological tissue regeneration. The TCC blocks the flexo-extension movement of the ankle during walking reducing the shear and torsion forces that constitute an additional factor that opposes tissue repair of the injury.

Effective off-load is the key to treating ulcers in diabetic foot syndrome, especially in case of neuropathic plantar lesions of the forefoot and midfoot. The sensory motor neuropathy in fact acts both by altering the anatomy of the foot creating areas of hyperload (metatarsal heads, mediatarsal region, apex of the fingers), and by altering the gait biomechanics (prolonged stationing in the metatarsal-phalanx regions).

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Along with assessment and management of infections and assessment of vascular condition, TCC is **an essential component for the treatment of diabetic ulcerative lesions** and represents the **most effective and reliable method** to download DFUs.

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Clinical evidence

The use of TCC as the most effective off-loading method is strongly supported by clinical evidence. Randomized clinical trials provide strong evidence of reduced healing time and greater numbers of healed DFUs for TCC when compared to other methods.

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Lavery et al. (2014) showed that TCC had a wound closure rate of 88.9%, compared to 50% for healing sandal and 40% for shear-reducing walking boot. They also showed that median time to healing for TCC was 5.4 weeks, compared to 8.9 weeks for a healing sandal and 6.7 weeks for sheer reducing walking boot.

Armstrong et al. (2001) determined that TCC had a healing rate of 88.9%, compared to 65% for a removable cast walker and 58.3% for a half-shoe, and time to healing was least for TCC at 4.8 weeks, compared to 7.2 weeks for the removable cast walker and 8.7 weeks for the half-shoe. This data was reiterated by Lavery almost 14 years later in a newer study (Figure 1).

Figure 1


Figure 1 - Randomized, controlled clinical trial from 2014, comparing total contact casting (TCC), healing sandal (HS) and shear-reducing walking boot (SRB). Data replotted from Lavery et al. 2014.[40]

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Contraindications

Total contact casting is contraindicated in the presence of untreated infection or osteomyelitis and in patients with severe peripheral arterial disease.

It should be used with caution in deep or heavily draining wounds and in ataxic, blind or severely obese patients.

Patient compliance is also necessary to minimize complications. Patients may require assistance with activities of daily living, such as bathing and sleeping. Moreover, any advanced wound-healing adjunctive therapies that require daily applications may not be suitable for use with patients using a non removable device, such as the TCC. In addition, using TCC has a moderate risk of causing new conflict lesions and is not indicated in claustrophobic patients.

*In situations requiring prolonged casting, **joint rigidity** and **muscular atrophy** have been documented. The patient should also be advised that certain designs of TCCs might exacerbate postural instability.*

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Charcot neuro-osteoarthropathy

*Charcot neuropathic osteoarthropathy (CN), as complication of diabetic neuropathy, is a disease which, **further to the progressive alteration of the bone structure of foot and/or ankle**, leads to the development of serious deformities with high risk of ulceration, infection and, as a consequence, to major amputation.*

*The **physiopathological process** may lead to progressive formation of deformities and subsequent ulcerations close to the pathological weightbearing areas. The development of infection and osteomyelitis, after the onset of an ulceration, may become a dangerous complication that leads to an increased risk of amputation unless an effective treatment, either non-surgical or surgical, is set in place. In literature ulceration is described in 40-50% of the patients affected by CN.*

*TCC is considered the **gold standard** not only to heal DFUs but also to protect the foot during the early vulnerable phases of Charcot disease when risk of bone fracture and dislocations is high.*

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An innovative structure for immobilization in total contact.

SCUDO IS A BRACE FOR IMMOBILIZATION OF THE TIBIO-TARSAL JOINT ABLE TO OBTAIN CONTROL OVER THE FOOT LOAD AND TOTAL CONTACT WITH THE PATIENT'S LIMB.

*FGP has revolutionized the constructive approach, conceiving the brace as a single piece, without removable parts, equipped with the innovative A.C.T., pneumatic system (P.I.S.), able to replicate the plantar conformation, **thus ensuring the immobilization, discharge and reduction of skin friction.***

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WHY CHOOSE SCUDO

- Once closed SCUDO offers a *scientifically comparable immobilization obtained with a cast.*
- Allows you to *customize the injury/wound monitoring path* based on the specific needs of the patient by simplifying the management of the brace.
- The *structure of SCUDO* makes it suitable for use also for those patients for whom it is contraindicated the use of a total contact cast.
- It *reduces the time and costs* necessary to replace the total contact cast that can be implemented only by specialized technical personnel.
- *Facilitates the operations of customization* on the specific needs of the patient by the technician.
- *Increases patient compliance* to treatment.
- The brace has been designed to be *managed independently* with methods and times indicated by the doctor.
- It offers *greater safety in walking* thanks to the shape of the sole and the non-slip material.

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TECNOCAST MAKES SCUDO A REVOLUTIONARY BRACE WITH TOTAL CONTACT TECHNOLOGY, PRESENTS CUTTING-EDGE CONSTRUCTION SOLUTIONS (INCLUDING THE ADVANCED CLAMPING TECH PATENT) THAT INTEGRATE WITH THE USE OF TECHNOLOGICAL MATERIALS TO PROPOSE A NEW APPROACH AMONG THE BRACES FOR THE DIABETIC FOOT.

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PATENTED LOCKING TENSION MODULATION SYSTEM ALLOWING PRECISION ADJUSTMENT BY:

- *rotating drive mechanism that requires less effort to obtain the correct closing tension;*
- *tension release mechanism through a simple "pull-out" movement of the knob;*
- *two different tension circuits for the upper and lower part that allow the brace to be adapted according to the specific needs of the patient.*

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PATENTED "CLAMP" CLOSURE SYSTEM, EXCLUSIVE FOR SCUDO, THAT ALLOWS TO OBTAIN THE BIOMECHANICAL ACTION OF LOCKING THE TIBIO-TARSICA JOINT THANKS TO:

- *patented double shell structure;*
- *high stiffness technopolymers;*
- *shell geometry designed in surface and in cuts that allows the brace to adapt to any conformation.*

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PNEUMATIC PLANTAR SPHERE SYSTEM TO OBTAIN TOTAL CONTACT WITH THE SURFACE OF THE FOOT BED BY:

- *circuit structure that uses the weight of the patient to assume the configuration of total contact with the foot and maintain it throughout the wearing period;*
- *possibility of restoring, every time you remove the brace, the neutral configuration, acting on the appropriate button of the circuit;*
- *presence of a structure that contains the insole that protects the patient from any overpressure.*

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- 1 HIGH STIFFNESS TECHNOPOLYMER STRUCTURE**
for features and shape ensures the stabilization and blockage of the tibio-tarsic joint and foot, thereby reducing friction.
- 2 LARGE MALLEOLAR OPENINGS**
with contoured edge surfaces that prevent risks of localized overpressure.
- 3 LARGE OVER HEEL OPENING**
with contoured edges to prevent localized overpressure.
- 4 SPIN LOCK SYSTEM**
to adjust the closing tension of the brace.
- 5 COUPLING SYSTEM BETWEEN PADDING AND RIGID STRUCTURE**
prevents the relative movement between the two elements once the brace worn.
- 6 NON-SLIP SOLE**
with boat shape; with angles designed to ensure proper discharge, stability and fluidity of movement.
- 7 LONGITUDINAL CUTS**
designed to ensure the necessary adaptability and stiffness.

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WHILE ITS DESIGN FOR FUNCTIONALITY WAS RIGOROUS, ATTENTION WAS PAID TO DESIGN OF THE PRODUCT SO THAT IT WOULD ALSO BE PLEASANT.

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Button for the opening of the circuit **PNEUSPHERE INSOLE SYSTEM**.

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PADDING in external gripping material and internal open cell material to ensure breathability, safety and comfort.

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TENSIONERS WITH CLAMPING CLOSURE facilitated by the presence of the pull and ergonomic shapes.

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*Laboratory test has shown that **SCUSO brace is comparable to total contact cast** in terms of plantar pressure distribution both in static and dynamic.*

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

TO CARRY OUT THE COMPARISON A SCUSO BRACE (FIGURE A) AND A TOTAL CONTACT CAST MADE BY SPECIALIZED PERSONNEL (FIGURE B) WAS USED.



ORTOSTASI / ORTHOSTATIC

GAMBALETTO GESSATO
TOTAL CONTACT CAST

Pressione
Pressure



TUTORE SCUDO
SCUDO BRACE

Area di contatto cm²
Contact area cm²



DEAMBULAZIONE / GAIT

GAMBALETTO GESSATO
TOTAL CONTACT CAST

Pressione
Pressure

TUTORE SCUDO
SCUDO BRACE

Area di contatto cm²
Contact area cm²

Misure (ambidestro) Measures (ambidextrous)

TAGLIE SIZE	CALZATA SCARPA SHOE SIZE
SMALL	35-38
MEDIUM	39-43
LARGE	44-46

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Indications

Brace designed to immobilize the foot and ankle.

Allows to maintain the surface of foot contact and limit skin stress. Useful, for example, in cases of:

- Acute phase treatment and secondary prevention of diabetic foot ulcers
- Treatment of Charcot's foot
- Treatment of vascular foot
- Treatment of fractures and amputations of the midfoot, forefoot and hindfoot

Technical characteristics

DIABETIC FOOT BRACE MADE WITH THE NEW TECNOCAST TECHNOLOGY BY FGP

STABILISING FRAME

- High rigidity technopolymer outer hull
- Patented double shell structure, ADVANCED CLAMPING TECH
- Non - slip sole with a boat-like shape and a predisposition to the positioning of the heel unloader wedges

CLOSURE SYSTEM

- Equipped with two independent patented SPIN LOCK systems for the closing tension modulation
- Operable via rotation for closing and pull - out for opening
- Equipped with tensioners that are easy to grip thanks to the ergonomic hook and quick release pull

PADDING

- Made by gripping material outside and open cell material inside for greater comfort and breathability of the patient
- With additional EVA pad (30 SHORE) applicable via gripping system in the brace pad calf area, marked by the petroleum blue seam (this optional element to be used only if volumetric calf compensation is required)

INSOLE

- Pneumatic plantar micro sphere system, PNEUSPHERE INSOLE SYSTEM, that's able to copy the patient's foot conformation for total foot contact
- The pneumatic circuit allows the correct configuration to be restored every time worn
- The insole is protected in a pocket with padded side protectors for a grater patient safety and comfort

NOTE

A set of wedges for the heel unload can be provided on request. The set consists of 2 wedges, 15 mm and 30 mm height, applicable to the sole with the dedicated screw set in the space marked on the sole by the hatching.

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Indications

Brace designed to immobilize the foot and ankle.
Allows to maintain the surface of foot contact and limit skin stress. Useful, for example, in cases of:

- Heel ulcer acute phase treatment
- Treatment of vascular foot
- Hind foot lesions post operative treatment

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Technical characteristics

DIABETIC FOOT BRACE MADE WITH THE
NEW TECNOCAST TECHNOLOGY BY FGP

STABILISING FRAME

- High rigidity technopolymer outer hull
- Patented double shell structure, ADVANCED CLAMPING TECH
- Non - slip sole with a boat-like shape
- Heel unloader wedge with height chosen by the technician (15 mm or 30 mm)

CLOSURE SYSTEM

- Equipped with two independent patented SPIN LOCK systems for the closing tension modulation
- Operable via rotation for closing and pull - out for opening
- Equipped with tensioners that are easy to grip thanks to the ergonomic hook and quick release pull

PADDING

- Made by gripping material outside and open cell material inside for greater comfort and breathability of the patient
- With additional EVA pad (30 SHORE) applicable via gripping system in the brace pad calf area, marked by the petroleum blue seam (this optional element to be used only if volumetric calf compensation is required)

INSOLE

- Pneumatic plantar micro sphere system, PNEUSPHERE INSOLE SYSTEM, that's able to copy the patient's foot conformation for total foot contact
- The pneumatic circuit allows the correct configuration to be restored every time worn
- The insole is protected in a pocket with padded side protectors for a grater patient safety and comfort

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Padding made by external gripping material and open cell material in the inside to ensure breathability, safety and comfort.

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0.75 mm stainless steel cable for the SPIN LOCK system.

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Additional pad made of EVA (30 Shore) applicable, by gripping system, to the brace padding for volumetric compensation of the calf.

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Set of wedges, composed of a wedge in height 15 mm and one in height 30 mm, made of technopolymer and attachable to the sole by means of special screws. Equipped with non-slip rubber at the base. Optionally usable on the brace with the aim of discharging the heel area.

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Complete replacement of the pneumatic system consisting of:

- microsphere insole
- cover
- pneumatic circuit

MISURE / MEASURES

S - M - L

CODICE D'ORDINE / ORDER CODE

DBT-RICPIS

Bibliografia / References

- 1 Armstrong, D.G., et al., It's not what you put on, but what you take off: techniques for debriding and off-loading the diabetic foot wound. *Clin Infect Dis*, 2004. 39 Suppl 2: p. S92-9.
- 2 Raspovic, A. and K.B. Landorf, A survey of offloading practices for diabetes-related plantar neuropathic foot ulcers. *J Foot Ankle Res*, 2014. 7: p. 35.
- 3 Snyder, R.J., et al., The management of diabetic foot ulcers through optimal off-loading building consensus guidelines and practical recommendations to improve outcomes. *J Am Podiatr Med Assoc*, 2014. 104(6): p. 555-67.
- 4 Armstrong, D.G., et al., Off-loading the diabetic foot wound: a randomized clinical trial. *Diabetes Care*, 2001. 24(6): p. 1019-22.
- 5 Lavery, L.A., et al., Reducing dynamic foot pressures in high-risk diabetic subjects with foot ulcerations: a comparison of treatments. *Diabetes Care*, 1996. 19: p. 818-821.
- 6 Lewis, J. and A. Lipp, Pressure-relieving interventions for treating diabetic foot ulcers. *Cochrane Database Syst Rev*, 2013. 1: p. Cd002302.
- 7 Schmidt, V.E., J.H. Somerset, and R.E. Porter, Mechanical properties of orthopedic plaster bandages. *J Biomech*, 1973. 6(2): p. 173-85.
- 8 Kowalski, K.L., J.D. Pitcher, Jr., and B. Bickley, Evaluation of fiberglass versus plaster of Paris for immobilization of fractures of the arm and leg. *Mil Med*, 2002. 167(8): p. 657-61.
- 9 Greenhagen, R.M. and D.K. Wukich, Total contact casting for neuropathic ulcers: a lost art? *J Diabetic Foot Complications*, 2009. 1(4): p. 2-9.
- 10 Fife C.E. et al., Diabetic foot ulcer off-loading: The gap between evidence and practice. Data from the US Wound Registry. *Adv Skin Wound Care*. 2014 Jul;27(7):310-6.
- 11 Capobianco, C.M., J.J. Stapleton, and T. Zgonis, The role of an extended medial column arthrodesis for Charcot midfoot neuroarthropathy. *Diabet Foot Ankle*, 2010. 1.
- 12 Wertsch, J.J., et al., Plantar pressures with total contact casting. *J Rehabil Res Dev*, 1995. 32(3): p. 205-9.
- 13 Armstrong, D.G., et al., Off-loading the diabetic foot wound: a randomized clinical trial. *Diabetes Care*, 2001. 24(6): p. 1019-22.
- 14 Morona, J.K., et al., Comparison of the clinical effectiveness of different off-loading devices for the treatment of neuropathic foot ulcers in patients with diabetes: a systematic review and meta-analysis. *Diabetes Metab Res Rev*, 2013. 29(3): p. 183-93.
- 15 Mueller, M.J., et al., Effect of Achilles tendon lengthening on neuropathic plantar ulcers. A randomized clinical trial. *J Bone Joint Surg Am*, 2003. 85-a(8): p. 1436-45.
- 16 Piaggese, A., et al., An off-the-shelf instant contact casting device for the management of diabetic foot ulcers: a randomized prospective trial versus traditional fiberglass cast. *Diabetes Care*, 2007. 30(3): p. 586-90.
- 17 Caravaggi, C., et al., Effectiveness and safety of a nonremovable fiberglass off-bearing cast versus a therapeutic shoe in the treatment of neuropathic foot ulcers: a randomized study. *Diabetes Care*, 2000. 23(12): p. 1746-51.
- 18 Lavery, L.A., et al., Randomised clinical trial to compare total contact casts, healing sandals and a shear-reducing removable boot to heal diabetic foot ulcers. *Int Wound J*, 2014.
- 19 Rogers L.C., Frykberg RG et al. The Charcot foot in diabetes. *Diabetes Care*. 2011; 34(9): 2123-9.

VISIT OUR WEBSITE
scudo.fgpsrl.it

FGP srl
Via Alessandro Volta 3
I-37062 Dossobuono VR
T +39 0458600867
F +39 0458600835
www.fgpsrl.it - fgp@fgpsrl.it

AZIENDA CERTIFICATA
ISO 13485:2016

Scarica l'App
FGPOrteSystem®