

SPARE PARTS



HULL PAD

Padding made by external gripping material and open cell material in the inside to ensure breathability, safety and comfort.

MEASURES

S - M - L

ORDER CODE

DBT-RICIMB



SPIN LOCK TENSIONING CABLE

0.75 mm stainless steel cable for the SPIN LOCK system.

MEASURES

S - M - L

ORDER CODE

DBT-RICCAVOSL



ADDITIONAL CALF PAD

Additional pad made of EVA (30 Shore) applicable, by gripping system, to the brace padding for volumetric compensation of the calf.

MEASURES

S - M - L

ORDER CODE

DBT-RICPAD



UNLOADING WEDGES

Set of wedges, composed of a wedge in height 15 mm and one in height 30 mm, made of technopolymer and attachable to the sole by means of special screws. Equipped with non-slip rubber at the base. Optionally usable on the brace with the aim of discharging the heel area.

MEASURES

S - M - L

ORDER CODE

DBT-RICCUNEI



COMPLETE PNEUSPHERE INSOLE SYSTEM

Complete replacement of the pneumatic system consisting of:

- microsphere insole
- cover
- pneumatic circuit

MEASURES

S - M - L

ORDER CODE

DBT-RICPIS