

SCUDO PLUS

TOTAL CONTACT DIABETIC CAST PLUS

ORDER CODE: DBT-SCUDOPLUS



Indications

Brace designed to immobilize the foot and ankle. Allows to maintain the surface of foot contact and limit skin stress. Useful, for example, in cases of:

- Heel ulcer acute phase treatment
- Treatment of vascular foot
- Hind foot lesions post operative treatment

Measures (ambidextrous)

SIZE	SHOE SIZE
SMALL	35-38
MEDIUM	39-43
LARGE	44-46

WEDGE
15 mm
30 mm

Technical characteristics

DIABETIC FOOT BRACE MADE WITH THE NEW TECNOCAST TECHNOLOGY BY FGP

STABILISING FRAME

- High rigidity technopolymer outer hull
- Patented double shell structure, ADVANCED CLAMPING TECH
- Non - slip sole with a boat-like shape
- Heel unloader wedge with height chosen by the technician (15 mm or 30 mm)

CLOSURE SYSTEM

- Equipped with two independent patented SPIN LOCK systems for the closing tension modulation
- Operable via rotation for closing and pull - out for opening
- Equipped with tensioners that are easy to grip thanks to the ergonomic hook and quick release pull

PADDING

- Made by gripping material outside and open cell material inside for greater comfort and breathability of the patient
- With additional EVA pad (30 SHORE) applicable via gripping system in the brace pad calf area, marked by the petroleum blue seam (this optional element to be used only if volumetric calf compensation is required)

INSOLE

- Pneumatic plantar micro sphere system, PNEUSPHERE INSOLE SYSTEM, that's able to copy the patient's foot conformation for total foot contact
- The pneumatic circuit allows the correct configuration to be restored every time worn
- The insole is protected in a pocket with padded side protectors for a grater patient safety and comfort